

EDGE 3

4.5 MPH iLEVEL® READY



Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. The Edge 3 is available in a total of 20 color options, including seven bold colors.



RASPBERRY BERET

Q
QUANTUM®
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic; 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 6.25 mph ¹¹ Up to 4.5 mph w/iLevel activated
GROUND CLEARANCE	2.7" (battery tray)
TURNING RADIUS ²	20.5"
OVERALL LENGTH ²	32.23" without front riggings
BASE WIDTH	24.4"
SEATING SIZES	
Synergy®/TRU-Balance® 3 Static Seat	Width: 10-22" Depth: 10-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Tilt and iLevel®	Width: 12-22" Depth: 12-22"
TRU-Balance® Flex Seating	Width: 15-22" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS ¹⁰	
Synergy®/TRU-Balance® 3 Static Seat ⁵	Static: 16.5 - 18.5" Synergy: 15.875 - 20.875"
TRU-Balance® 3 Power Tilt ⁵	17.5 - 19.5"
TRU-Balance® 3 Tilt & iLevel® ⁵	17.5 - 19.5"
TRU-Balance® Flex Seating	16.5 - 18.25"
MANUAL RECLINE ⁶	Yes
BATTERY SIZE ³	NF-22 (2)
BATTERY WEIGHT ⁷	38 lbs. each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel ⁶
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	138.6 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856/K0857), Multiple Power (K0861)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability
- 10) 16.5" seat-to-floor height is not available with power tilt or iLevel®. 19.5" seat-to-floor height is not available with static seat or power recline.
- 11) 6.25 mph speed available only with Q-Logic 3.

Features

- Upgraded SRS (Smooth Ride Suspension)
- 4.5 mph with iLevel®
- Fender lights and USB charger included
- Bold new colors
- Made in America of foreign and domestic parts

7 Bold Colors*



BLUE
JEAN

CHERRY
BOMB

LEMON
CRUSH

GREEN
MACHINE



RASPBERRY
BERET



EVOLUTION
ORANGE



SUGAR
PLUM

HYDROGRAPHICS



CARBON
FIBER

DIGITAL
CAMO

FADED
GLORY

* Visit quantumrehab.com to see all 20 available colors.



PROUDLY MADE IN
AMERICA 
OF FOREIGN AND DOMESTIC PARTS

QUANTUM®

#1 FOR REHAB POWER

(US) 866-800-2002

(Canada) 888-570-1113

quantumrehab.com

