Features, Performance, and Travel Ready






The Go-Go® Sport delivers high-performance and convenience on the go. With a 325 lb. weight capacity, a charger port in the tiller, standard front LED lighting and a maximum speed up to 4.7 mph., the Go-Go® Sport is feature rich and travel ready.



Wraparound Delta Tiller







Specifications

3W: SC73 (K0801) 4W:SC74

Weight Capacity	325 lbs.
Maximum Speed ¹	Up to 4.7 mph
Ground Clearance ²	1.5" at motor
Turning Radius	3W: 36″ 4W: 50.5″
Overall Length ²	3W: 42″ 4W: 43″
Overall Width ²	21.25″
Front Tires	3W: 3" x 8" solid 4W: 3" x 9" solid
Rear Tires	3″ x 9″ solid
Range Per Charge ^{1,4}	3W: Up to 13 miles 4W: Up to 12 miles
Total Weight Without Batteries	3W: 93 lbs. 4W: 99 lbs.
Heaviest Piece When Disassembled	3W: 35 lbs. (front section) 4W: 44 lbs. (front section)
Standard Seating	<i>Type:</i> Foldable, compact, vinyl-covered, molded plastic Weight: 23 lbs. Material: Black Vinyl Dimensions: Width: 17" Depth: 17"
Motor Output	24-volt DC motor, sealed transaxle
Brakes	Regenerative and electromechanical
Battery Requirements ^{3,5}	(2) 12V deep-cycle (standard) Size: (2) 18 AH, 14 lbs. each Weight: 29.5 lbs. per battery pack
Battery Charger	Off-board, 2 amp
Colors	2 sets with every scooter!
	Red Blue

- Charger port on tiller
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Deck basket for additional storage (4-wheel only)
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Front tiller-mounted basket
- Includes two sets of easily changeable, red and blue colored panels
- Long-lasting, high-intensity, front LED headlight offers optimal pathway illumination
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tires

Accessories

- Additional 18AH battery box
- Cane/crutch holder
- Cell phone holder
 - Cup holder
- 50" Lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger



- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 2 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).
- AGM or gel cell type required.

| (CAN) 888-570-1113 | pridemobility.com | f 🖸 💟 🞯

- Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5 Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or -3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



Pride offers a variety of lifts to **transport your mobility product with ease.** To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at <u>pridemobility.com/lift-finder/</u>

(US) 800-800-8586