



3 wheel GB118 4 wheel GB148

Seat Color Black Vinyl

Body Colors Red BI

The all new Golden BuzzAround Extreme Scooter features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full size scooter! The BuzzAround Extreme Scooters come standard with 2-U1 batteries - industry leading and the largest battery in this category - for maximum operating range up to 18 miles. The Extreme also features a brand new ultrabright LED headlight and taillight, 330 lb. weight capacity, larger 9" tires, and comfort suspension springs for the smoothest ride. Go the extra mile on a BuzzAround Extreme Scooter!





Item	GB118	GB148
Description	EX	EX
Medicare Code	Not Coded	Not Coded
Weight Capacity	330 lbs.	330 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	5 mph	5 mph
Operating Range up to <sup>1</sup>	18 miles	18 miles
Jnder Deck Clearance	4"	4"
Turning Radius	42 "	56"
Type Batteries	2-U1	2-U1
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brake	es Yes	Yes
Charger	Off-board	Off-board
Controller	90 Amp Dynamic R-Series	90 Amp Dynamic R-Series
Dimensions		
Length	45"	44"
Width	22"	22"
Height (ground to top of back of seat)	35"-37"	35"-37"
Ground to Top of Seat	22.75"-24.75"	22.75"-24.75"
Ground to Top of Deck	6.75"	6.75"
Deck to Top of Seat	16"-18"	16"-18"
Number of Seat Height P	ositions 3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	34.7"	33.8"
Weight of Unit: (assembled)	148 lbs.	155 lbs.
Front Half	46 lbs.	53 lbs.
Rear Half	30 lbs.	30 lbs.
Seat w/Arms	23 lbs.	23 lbs.
Batteries	24 lbs. per pack	24 lbs. per pack
Tires:	Solid	Solid
Front	9"	9"
Rear	9"	9"
Rear Anti-Tip Wheels	2"	2"
Standard Seat	Stadium Style	Stadium Style
Seat Back Height - no headrest	16"	16"
Width x Depth	18" x 16"	18" x 16"
Color	Black Vinyl	Black Vinyl
Optional Seat:	Stadium Style	Stadium Style
Back Height - no headrest	16"	16"
Width x Depth	20"x 18"	20"x18"
Color	Black Vinyl	Black Vinyl
Warranty:	2.2	
Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics <sup>2</sup>	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

<sup>&</sup>lt;sup>1.</sup> Battery range tested at 200 lbs. and will vary due to rider weight, drive surface and terrain.

## Wireless disassembly in 4 easy steps!



**Step 1:** Pop off the seat!



Step 2: Pull out the battery pack!



**Step 3:** Fold down tiller and lock in place!



**Step 4:** Pull forward on the drivetrain release lever to separate the front and rear frame sections!



Disassembles into 6 pieces!

4 wheel model shown


## **Features**

- 330 lb. weight capacity
- 4" ground clearance
- Full size & portable
- Larger 90 amp controller
- Larger 320 watt motor
- Larger 9" tires
- U1 battery for maximum range
- Ultrabright LED tail light

- Comfort suspension springs
- Full size, yet portable & easy to disassemble



A battery charging port is conveniently located on the tiller under the control panel.





<sup>&</sup>lt;sup>2.</sup> Electronics warranty excludes batteries.