

Invacare® Pronto® M61™ Power Wheelchair

POWER  
CUSTOM  
REHAB  
SEATING  
SUPPORTS

# Pronto® M61™



Yes, you can.®





Yes, you can.®

Versatility comes standard with the Invacare® Pronto® M61™ power wheelchair with SureStep® and power elevating seat. The M61 boasts a power elevating seat that provides up to five inches of elevation with six wheels on the ground for optimum stability. In the elevated position, the M61 provides smooth driving and accessibility.\* The optional captain's seat offers a plush charcoal grey vinyl semi-recline seat and adjustable height arms, providing optimal consumer comfort. The M61-R comes standard with Rehab (ASBA) seating. The Pronto M61 with SureStep provides a smooth ride over transitions and thresholds up to two inches while maintaining stability. The M61 design offers sleek styling and comfort while optimizing maneuverability and accessibility. \*Do not drive on an incline in the elevated position.

**Invacare Corporation**  
www.invacare.com

## USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

## Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

© 2008 Invacare Corporation. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

Form No. 04-201  
Rev. 09/08



Shown with charcoal grey semi-recline seat.

## Pronto M61 with Van Seat

- Power elevating seat with up to five inches of elevation
- SureStep technology allows for smooth driving over thresholds and transitions
- Stylish design with choice of Candy Red or Deep Blue shroud colors
- Compact size for the tightest overall turning radius
- True center-wheel drive system for intuitive driving
- Six wheels on the ground provide optimum stability
- Fully programmable MK6i™ electronics for maximum function and precise driver control
- Adjustable joystick mount
- Light grey medium back van seat with flip-up adjustable-width arms
- Rehab (ASBA) seat available on M61-R
- Charcoal grey semi-recline seat (optional) with standard flip-up adjustable width arms
- Powerful 350-watt, high-torque motors for reliable performance
- Regenerative braking and automatic motor locks for dependable stopping
- 3-amp on-board charger for convenience
- 300 lb. weight capacity

### Center-Wheel Drive



Center-wheel-drive offers the tightest turning radius, as well as improved tracking and ease of driving, over all other Invacare configurations.

### SureStep®



This patented automatic system ensures that everyday obstacles and thresholds are handled elegantly for a smooth ride.

### Optional Rehab Seating Available (ASBA) on M61-R



Customize the M61-R for consumers with moderate seating and positioning needs.

## Specifications (w/standard light grey med. back seat)

Overall dimensions	
Length (with footboard folded)	34"
Length (with footboard down)	38" - 40"
Width (excluding joystick)	24"
Height	35" - 37"
Power seat elevation	5"
Seat dimensions	
Width	16", 18", 20"
Depth	16" - 20"
Weight	
With batteries	202 lb.
Without batteries	155 lb.
Wheel and caster sizes	
Front casters	6" x 2"
Drive wheels (flat-free)	10.5" x 3.5"
Rear casters	6" x 2"
Maximum speed	
	4.25 mph
Incline capability	
	9°
Weight capacity	
	300 lb.
Turning radius	
Front (with footboard down)	19.5"
Rear	19.5"
Ground clearance	2.6"
Limited warranties	
Base frame	60 months
Seat frame	36 months
Electronics	12 months
Motor/gearbox	12 months