

The Jazzy® 600 ES features patented Mid-Wheel 6® technology and Active-Trac® ATX Suspension for an extraordinary level of performance, maneuverability and the tightest turning radius in the Jazzy line. The Jazzy 600 ES offers true all-around performance indoors and out making it an excellent power chair choice.

#### **LARGE 14" DRIVE WHEELS**





SPHERICAL OMNI-CASTERS



ACTIVE-TRAC® ATX SUSPENSION

# AMERICA'S FAVORITE POWER CHAIR

Learn more online.









go**jazzy**.com

## Specifications

Weight Capacity 300 lbs.  Maximum Speed Up to 4 mph  Ground Clearance 2.75"  Turning Radius 20.5"  Overall Length 35.5"  Overall Width 24.5"  Drive Wheels 14"  Front Wheels 6" solid OMNI-casters
Ground Clearance 2.75"  Turning Radius 20.5"  Overall Length 35.5"  Overall Width 24.5"  Drive Wheels 14"
Turning Radius 20.5"  Overall Length 35.5"  Overall Width 24.5"  Drive Wheels 14"
Overall Length 35.5"  Overall Width 24.5"  Drive Wheels 14"
Overall Width 24.5"  Drive Wheels 14"
Drive Wheels 14"
Front Wheels 6" solid OMNI-casters
Rear Wheels 6" solid OMNI-casters
Suspension Active-Trac® (ATX)
<b>Drivetrain</b> Two-motor, in-line, mid-wheel drive
Braking System Regenerative and electro-mechanical
Standard Electronics 50A, Dynamic Shark controller
Battery Charger Off-board, 5 amp
Per-Charge Range Up to 17 miles
Battery Requirements (2) 12 volt, deep cycle, U-1
Battery Weight 25 lbs. each
Base Weight 124 lbs.
<b>Standard Seat Weight</b> 52 lbs. (high-back)
<b>Standard Seat</b> 18" x 18-20" and 20" x 18-20"
Max. Seat Dimensions 20" x 20"
Warranty Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## Standard Color



Cardinal Red

### **Features**

- · Depth adjustable, black high-back seat
- · Universal, fixed, in-line joystick mount
- · Flip-up height adjustable foot platform
- · Black, flat-free, non-marking tires
- 60" lap belt

## Options

- · Synergy® Seat
- · Solid seat pan
- Cup holder
- · Weather cover
- Rear basket
- Oxygen tank holder
- · Walker holder
- Cane/crutch holder
- · Saddle bag

- · Heel loops
- Lap belt (50", 70", 80", 90")
- · Swing-away footrest
- · Elevating leg rests
- Military patches
- · Swing-away joystick
- · XLR USB charger

## Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Zeus
- Commander
- Outlander Full Platform



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder

## Models

Jazzy 600 ES 2S-SS Code: K0822 Jazzy 600 ES 2S-C Code: K0823











