

# Jazzy<sup>®</sup>

## SELECT 6

The Jazzy Select<sup>®</sup> 6 incorporates sporty design and performance features for a truly unique ride experience. Pride's patented Active-Trac<sup>®</sup> Suspension combined with six wheels on the ground and an innovative articulating beam allow for superior performance indoors and out, with added stability over uneven terrain and ramp transitions.



### ACTIVE-TRAC<sup>®</sup> SUSPENSION



**JAZZY ARMOR<sup>®</sup> AND  
COLOR-THROUGH  
SHROUD**



**COMFORT HIGH-BACK SEAT  
WITH REPLACEABLE BACK  
AND SEAT COVERS**

AMERICA'S FAVORITE  
POWER CHAIR

Learn more online.



[gojazzy.com](http://gojazzy.com)

## Specifications

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed</b>	Up to 4 mph
<b>Ground Clearance</b>	1.75"
<b>Turning Radius</b>	21.75"
<b>Overall Length</b>	35.5"
<b>Overall Width</b>	23.5"
<b>Drive Wheels</b>	10"
<b>Front Wheels</b>	6" solid casters
<b>Rear Wheels</b>	6" solid casters
<b>Suspension</b>	Active-Trac® (ATX)
<b>Drivetrain</b>	Two-motor, in-line, mid-wheel drive
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Standard Electronics</b>	40 amp, PG GC 3 controller
<b>Battery Charger</b>	Off-board, 3 amp
<b>Per-Charge Range</b>	Up to 15 miles
<b>Battery Requirements</b>	(2) 12 volt, deep cycle, U-1
<b>Battery Weight</b>	25 lbs. each
<b>Base Weight</b>	100 lbs.
<b>Standard Seat Weight</b>	43 lbs. (comfort, high-back)
<b>Standard Seat</b>	18" x 18-20"
<b>Max. Seat Dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## Colors



Jazzy Red



Jazzy Blue

## Features

- Jazzy® Armor® includes shroud and controller guards to protect against daily wear and tear
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## Options

- Synergy® Seat
- Solid seat pan
- Power elevating seat (non-swivel)
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick
- XLR USB charger

## Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

[pridemobility.com/Lift-Finder](http://pridemobility.com/Lift-Finder)

## Models

<b>Jazzy Select 6-SS</b>	Code: K0822
<b>Jazzy Select 6-C</b>	Code: K0823