

Jazzy[®]

SELECT ELITE

The Jazzy Select[®] Elite delivers high performance in a compact, stylish design for your active lifestyle. In-line, front-wheel drive technology gives the Jazzy Select Elite excellent stability and provides a tight turning radius around corners for optimal indoor maneuverability. The Jazzy Select Elite comes standard with Pride's exclusive Jazzy[®] Armor[®] for today's rental environment.



ADJUSTABLE FOOT PLATFORM



**JAZZY ARMOR[®] AND
COLOR-THROUGH
SHROUD**



**COMFORT HIGH-BACK SEAT
WITH REPLACEABLE BACK
AND SEAT COVERS**

**AMERICA'S FAVORITE
POWER CHAIR**

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.5"
Turning Radius	24.75"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 15 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	25 lbs. each
Base Weight	69 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Jazzy Red



Jazzy Blue

Features

- Color-through shroud
- Jazzy® Armor® includes shroud and controller guards to protect from daily wear and tear
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Options

- Synergy® Seat
- Solid seat pan
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Heel loops
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Cup holder
- Military patches
- Swing-away joystick
- XLR USB charger

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder

Models

Jazzy Select Elite-2S-SS

Code: K0822

Jazzy Select Elite-2S-C

Code: K0823