



3-Wheel shown with optional 22" wide high-back seat

The Victory® 10 series offers added style, comfort and performance for your active lifestyle. Standard features include LED pathway light, backlit battery gauge, feather touch disassembly, wraparound delta tiller, 10" solid tires and a front basket. A 400 lb. weight capacity, a per charge range of up to 15.5 miles and a maximum speed up to 5.25 mph compliment the versatility of the Victory 10. Available in Viper Blue and Candy Apple Red.













Specifications 3W:SC610 4W:SC710

| Weight Capacity | 400 lbs. |
|--|--|
| Maximum Speed ¹ | Up to 5.3 mph |
| Ground Clearance ³ | 3" at mid-frame |
| Turning Radius ³ | 3W: 42" 4W: 54" |
| Overall Length ³ | 3W: 45.25" 4W: 47" |
| Overall Width ³ | 22.25" |
| Seat-to-Ground Height | 20.5"-22" |
| Seat-to-Deck Height | 17.5″-19″ |
| Front Tires | 3.5" x 10" solid |
| Rear Tires | 3.5" x 10" solid |
| Suspension | N/A |
| Range Per Charge ^{1,2} | Up to 16 miles |
| Total Weight Without Batteries ^{3,7} | 3W: 133 lbs. 4W: 138 lbs. |
| Heaviest Piece When Disassembled ³ | 3W: 54 lbs. (front section) 4W: 61 lbs. (front section) |
| Standard Seating | Foldable molded plastic Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" width, 17" depth |
| Drive System | Rear-wheel drive sealed transaxle with 24-volt DC motor |
| Dual Braking System | Electronic, regenerative and electromechanical |
| Battery Requirements ^{3,4,5} | Size: (2) U-1(12V deep-cycle) standard, (2) 40 AH (optional) Weight: (U-1) 24.5 lbs. each, (40 AH) 32.5 lbs. each |
| Battery Charger | 3.5-amp, off-board (U-1), 5-amp, off-board (40 AH) |
| Colors | Candy Apple Red Viper Blue |

Features

- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Delta tiller with wraparound handles
- Backlit battery gauge
- Exclusive stylish, lightweight, non-scuffing, black, low-profile tires
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight pinchless hinge seat eight position rotation, sliders and viscoelastic foam inserts for increased comfort
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port
- Auto-connecting front-to-rear harness

Accessories

- Forearm crutch holder
- Cane/crutch holder (single, double)
- Oxygen tank holder
- Cup holder
- Rear basket

- Safety flag
- Walker holder
- Weather cover
- Saddlebag
- Cell phone holder
- XLR USB charger

ptions

40 AH batteries

ptional Seating

- 18"W x 18"D black vinyl reclining high-back seat, with sliders, height adjustable headrest, flip-up, height and width adjustable armrest
- 22"W x 18"D black vinyl reclining high-back seat height adjustable headrest, flip-up, height and width adjustable armrest
- Optional seat sizes available

Military Patches





Pride offers a variety of lifts to transport your mobility product with ease. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at pridemobility.com/lift-finder/

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery
- condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 AGM or gel cell type required.
- 5 Battery weight may vary based on battery manufacturer.
- 6 Includes standard seat weight.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.









