

# VICTORY<sup>®</sup> SPORT

4-WHEEL SCOOTER

## AN ACTIVE LIFESTYLE DESERVES A SPORTY RIDE

The **Victory<sup>®</sup> Sport** incorporates a sporty design with extra comfort and performance features for a truly unique ride experience. With a maximum speed of 8 mph — the fastest within the Victory Series — this scooter will accelerate your active lifestyle. Standard features include feather touch disassembly, wraparound delta tiller, 10" solid tires, front and rear suspension, high-back reclining seat, directional signals, a high-intensity LED lighting package, backlit battery gauge and a front basket. Available in Viper Blue and Candy Apple Red.



Candy Apple Red



Viper Blue



[facebook.com/PrideMobilityProductsCorp](https://facebook.com/PrideMobilityProductsCorp)





[youtube.com/PrideMobility](https://youtube.com/PrideMobility)

PRIDE MOBILITY PRODUCTS CORP. 10000 W. 10TH AVENUE, DENVER, CO 80202  
**PRIDE<sup>®</sup> MOBILITY SCOOTERS**  
**LIVE YOUR BEST<sup>®</sup>**

[pridemobility.com](https://pridemobility.com)

# SPECIFICATIONS

## VICTORY® SPORT 4 Wheel: SC710 DXW

Weight capacity:	350 lbs.
Maximum speed <sup>1</sup> :	Up to 8 mph
Ground clearance <sup>3</sup> :	2.25"
Turning radius:	51.75"
Overall length <sup>3</sup> :	47"
Overall width <sup>3</sup> :	26"
Seat-to-deck height <sup>3</sup> :	17.25" - 18.75"
Seat-to-ground height range:	21.5" - 24.25"
Front tires:	3.5" x 10" solid
Rear tires:	4" x 10" solid
Range per charge <sup>1,2</sup> :	Up to 11 miles
Total weight without batteries <sup>1,3</sup> :	173 lbs.
Heaviest piece when disassembled:	68 lbs. (front section)
Standard seating:	Type: Reclining high-back with sliders, height adjustable headrest, flip-up height and width adjustable armrest Weight: 45.5 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Motor output:	Rear-wheel drive, sealed transaxle with 24-volt DC motor
Brakes:	Electronic, regenerative and electromechanical
Battery requirements <sup>4</sup> :	(2) 12V deep-cycle (standard) Size: 40 AH Weight <sup>5</sup> : 32.5 lbs. each
Battery charger:	Off-board, 5 amp
Color:	  Candy Apple Red    Viper Blue

## FEATURES

- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Exclusive stylish, lightweight, non-scuffing, black, low-profile wheels
- Backlit battery gauge
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Black vinyl reclining high-back seat with sliders, height adjustable head rest, flip-up, height and width adjustable armrest
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port

## OPTIONS

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag
- Rearview mirror
- Lap belt
- Cell phone holder
- Touch-up paint
- Optional 22"x18" seat

## Pride® Lifts & Ramps

### Compatible Lifts:

- Backpacker® Plus
- Backpacker® AVP
- Olympian
- Outlander Full Platform
- Outlander DE
- Outlander XL
- Commander 400

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.  
<sup>2</sup> Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.  
<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.  
<sup>4</sup> AGM or gel cell type required.  
<sup>5</sup> Battery weight may vary based on battery manufacturer.  
<sup>6</sup> Includes standard seat weight.

### NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

VICTORYSPORT-SS-02-10-16



PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113

[pridemobility.com](http://pridemobility.com)